



Wahl® Powerful Heated Massager Guide

Attention: If the unit stops running as a safety precaution, please allow it to cool before attempting to use.



With the switch in the OFF position, plug the power cord into an AC 120-V outlet.



Place desired attachment head on stem of massager. Heads easily press on and pull off.



If heat or cold is desired, the special attachments must be used.

Note: See "Special Instructions"



Push the top switch for the heat setting. Use heat only with the Deep Heat attachment.



Begin at lowest speed setting. Gradually increase to desired massage speed by rolling dial toward top of massager. Roll dial toward bottom of massager to decrease speed.



Apply the massager to the targeted area, using continual movement. A circular motion is often most effective.

NOTE: See "Massage Techniques"

Wahl® Powerful Heated Massager Key Features

MASSAGE TECHNIQUES

SHOULDERS AND NECK: Massage each shoulder toward back of neck to point just behind and below ear. Then back to shoulder. Repeat several times.

ARMS: Start with a gentle massage of the palms and fingers. Next, using a circular motion, move upward to forearms. Then up to biceps and shoulders.

LEGS: If desired, massage bottoms of feet. Then, using a circular motion, massage each leg from calf on up and over thigh.

CHEST: Massage from center of chest outward, and up toward each shoulder.


BACK: Massage from top of spine to bottom, and from center outward.

ABDOMEN: Relax abdominal muscles by massaging from center outward, using a gentle, circular motion.

JOINTS: Relieve tired, aching knees and elbows with smooth, circular strokes.

FACE: Gently massage from the nose outward, taking care to avoid contact with the eye.

SCALP: Start at the back of the neck and massage up to the top of the head. Using this up and down motion, massage across the scalp.



****Colors and decoration styles may vary. Some kits may not include all attachments.****

If you are not completely satisfied with your Wahl Massager, please contact us before returning this product at WahlUSA.com or **1-800-767-9245**.