IMPORTANT SAFEGUARDS

When using your electrical appliance, basic precautions should always be followed, including the following:

READ ALL INSTRUCTIONS AND SAFEGUARDS BEFORE USING THIS MASSAGER.

DANGER

To reduce the risk of death or injury by electric shock:

1. Do not reach for an appliance that has fallen into water. Unplug immediately.
2. Do not use while bathing or in the shower.
3. Do not place or store appliance where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid.
4. Always unplug this appliance from the electrical outlet immediately after using and before cleaning.
5. Never use pins or other metallic fasteners with this appliance.
6. Carefully examine the covering before each use. Discard the appliance if the covering shows any sign of deterioration, such as checking, blistering, or cracking.
7. Keep Dry – Do not operate in a wet or moist condition.

WARNING

To reduce the risk of burns, electrical shock or injury to persons:

1. This appliance should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
2. Do not operate this appliance under a blanket or pillow. Excessive heating may occur and cause fire, electric shock or injury to persons.
3. Close supervision is necessary when this appliance is used by, on, or near children or disabled persons.
4. This massager contains a motor. To avoid risk of injury, it must not be used within six inches of an implantable cardiac device.
5. Use this appliance only for its intended use, as described in this manual. Do not use attachments not recommended by the manufacturer.
6. Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return to Wahl for examination and repair.
7. Do not carry this appliance by supply cord or use cord as a handle.
8. Keep the cord away from heated surfaces.
9. Never drop or insert any object into any opening on this appliance.
10. Do not use outdoors.
11. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
12. To disconnect, turn all controls to the off position, then remove plug from outlet.
13. Do not use massager in close proximity to loose clothing or jewelry.
14. Keep long hair away from massager while in use.
15. Consult your physician before using this massager in case of pregnancy, varicose veins, recent surgery, unexplained pain, or illnesses including but not limited to heart condition, diabetes, or osteoporosis.
16. Do not use massager on:
 - a. Sensitive skin, swollen or inflamed areas, or skin eruptions;
 - b. Areas of poor circulation, or insensitive skin;
 - c. The front of the neck;
 - d. Genital areas;
 - e. Sleeping or unconscious persons.
17. Do not use this massager for more than 20 minutes in any one area.
18. Before you unplug from outlet, hold the power button for 3 seconds to turn OFF and press the heat button to turn OFF, then remove plug from outlet.
19. Temperatures sufficiently high to cause burns may occur regardless of the control setting. Do not use on an infant or disabled person or on a sleeping or unconscious person. Do not use on insensitive skin or on a person with poor blood circulation. Check the skin in contact with the heated area of the appliance frequently to reduce the risk of blistering.
20. When storing this product, do not crush and avoid sharp folds.
21. Never operate the appliance with the air openings blocked. Keep the air openings (if any) free from lint, hair or the like.
22. This appliance is not intended for use by persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Keep children under the age of 13 away from this machine.

Some of our massagers have heat attachments and gel packs.

For a unit that has heat:

23. Use heated surfaces carefully. May cause serious burns. Do not use over insensitive skin areas or in the presence of poor circulation. The unattended use of heat by children or incapacitated persons may be dangerous.

For a unit that has a gel pack:

24. Never apply the frozen gel pack directly to skin as direct application may cause frostbite. Never apply the heated gel pack directly to skin as direct application may cause burns.

SOME UNITS MAY HAVE POLARIZED PLUGS

To reduce the risk of electric shock, this appliance has a polarized plug (one blade is wider than the other). This plug will fit in a polarized outlet only one way. If the plug does not fit fully into the outlet, reverse the plug. If it still does not fit, contact a qualified electrician to install the proper outlet. Do not change the plug in any way.

CLEANING/STORAGE

Always unplug massager before cleaning.

To clean, simply wipe unit with a dry, soft cloth. Do not immerse in water or use any harsh cleaners or solvents. When storing unit, make sure cord is not twisted or kinked.

SERVICE AND REPAIR

This unit has no user-serviceable parts. If you have problems with your massager, **do not return it to store where purchased.** Send it prepaid with a note explaining the problem to:

**Wahl Clipper Corporation
Attention: Repair Department
3001 Locust Street
Sterling, Illinois 61081
Telephone: (800) 767-9245**

For products purchased outside the United States, contact your Wahl distributor.

1-YEAR LIMITED WARRANTY

VALID IN U.S.A. ONLY

If your Wahl product fails to operate satisfactorily within one (1) year from the date of original purchase or receipt as a gift, Wahl will, at our option, repair or replace the product and return it to you **FREE OF CHARGE. Do not return this product to store where purchased.** Products needing repair under this warranty should be returned with postage pre-paid to the following address or call 1-800-767-9245 for assistance: **Wahl Clipper Corporation, Attn: Repair Department, 3001 Locust Street, Sterling, Illinois 61081.**

When mailing, please include a note describing the defect and the date of original purchase or receipt as a gift. Products out of warranty will be repaired at our standard repair charges.

This warranty does not cover products damaged by:

1. Modification or repair by anyone other than Wahl Clipper Corporation or a Wahl authorized service center;
2. Commercial (non-consumer) use or unreasonable use.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state. No other written express warranty is given. No responsibility is assumed for incidental or consequential damages of any type or for the use of unauthorized attachments. Some states may not allow the exclusion or limitation of incidental or consequential damages, so the limitation may not apply to you. Warranty valid in U.S.A. only. Outside U.S.A., check with store for Wahl authorized distributor's warranty.

REIMBURSEMENT

This unit may qualify for reimbursement through your work FSA (flexible spending account), Medicare, Medicaid, and some insurance companies.

For Products Marked with an FCC Approval:

1. This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
2. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this device.
3. NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
 - Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help.

**SAVE THESE INSTRUCTIONS
FOR HOUSEHOLD USE ONLY**

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